

GYPB-glycophorin B (MNS blood group) Gene

Reference: PTXBC069310

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GYPB/GYPA fusion; GYPB-A-B fusion; GYPB-A fusion; CD235b; GPB; GYP; MNS; PAS-3; glycophorin-B; SS-active sialoglycoprotein; Ss blood group; blood group system MNSs, Mur-like; blood group system MNSs, St(a) type E; blood group system MNSs, St(a) type F; gl		

For research use only.

GENE DESCRIPTION	glycophorin B (MNS blood group)
GENE ID	2994
ORIGIN SPECIES	Human
ACCESSIONS	BC069310
VECTOR	pDONR223
NCBI SYMBOL	GYPB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatggaaaaataatctttgtattactattgtcagaaattgtgagcatatcagcattaaagtaccactgaggtggcaatgcacacttcaacctcttctcagtcacaaagagttacatctcatcacagacaaatggagaaacgggacaaactgtccatctgttccactgtaccagctcctgtagtata
atactcattattttgtgtgatggctggtattattggaacgatcctcttaatttctacactattcgcgactgataaaggcatga