

IL20-interleukin 20 Gene

Reference: PTXBC069311

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-20; IL10D; ZCYTO10; interleukin-20; cytokine Zcyto10; four alpha helix cytokine; interleukin 20		

For research use only.

GENE DESCRIPTION	interleukin 20
GENE ID	50604
ORIGIN SPECIES	Human
ACCESSIONS	BC069311
VECTOR	pDONR223
NCBI SYMBOL	IL20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagcctctagcttgccttcagccttctctctgctgctgtttatctctatggactcctccactggactgaagacactcaattgggaagctgtgtgatccacaaacctcaggaaatcgaaatggattttctgagataggggcagtggtgcaagccaaagatggaaacattgacatca
gaatcttaaggagagactgagctcttgcagacacaaagcctgcgaatcgatgctgctcctcgcgcaatttgctaaagactctatctggacagggtatttaaaactaccagacctgacattatactctccggaagatcagcagcctgccaattcttcttaccatcaagaaggacctcgg
ctctgtcatgccacatgacatgccattgtggggaggaagcaatgaagaatacagccagattctgagtcatttgaagagctggaacctcaggcagcagttgtgaaggcttgggggaactagacattcttgcattggaaggagagacagaatag