

VGLL2-vestigial like 2 (Drosophila) Gene

Reference: PTXBC069316

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VGL2; VITO1; transcription cofactor vestigial-like protein 2; Vestigial and Tondu related protein 1; vestigial like 2; vestigial like family member 2		

For research use only.

GENE DESCRIPTION	vestigial like 2 (Drosophila)
GENE ID	245806
ORIGIN SPECIES	Human
ACCESSIONS	BC069316
VECTOR	pDONR223
NCBI SYMBOL	VGLL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctgtctggatgtatgtaccaagtctatgtctccgcagccctacttgcagccgcctacacccctaccaccagaaactagcctatttccaaatgcaggagcgcaggagtgcaatgccagccccagcagcagtgccagcggcagctctctatttccagccaaacccagccagtataaaagaggagaaggcagcccagagaagagcggccaccagaggcagagctacatcaactccgctgctctcttcttcttatttccaggggacatcagctccgtgggtggatgaacatttcagcaggccctgagccaacccagcagctactctctagctgtaccagcagaagcaccaaggagctctggccctggcgagctctgtcttattccctctgtgggtcatccctctgagctga