

MS4A6E-membrane-spanning 4-domains, subfamily A, member 6E Gene

Reference: PTXBC069322

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	membrane-spanning 4-domains subfamily A member 6E; membrane-spanning 4-domains, subfamily A, member 6E; membrane spanning 4-domains A6E		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 6E
GENE ID	245802
ORIGIN SPECIES	Human
ACCESSIONS	BC069322
VECTOR	pDONR223
NCBI SYMBOL	MS4A6E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacatcacacattttccaatgagaccatcataatgctcccatcaaatgtcatcaacttctcccaagcagagaaacccgaaccaccaaccagggcaggatagcctgaagaaacgtctacaggcaaaagtcaaagtattgggtgcatagcagcctggctggaagcattctgagtgc
tctgtctgccctgggtggtttcattctctgtctgtcaacccggctgcattaaatctgcctcattgcagtgtaagttggacgaaaaggatataccaaccagacttcttcttattgattatcattcaccttacacatggactgcatagagccaaagccagtctggctggaactctgtctctgatgc
tggtttctactgtgtggagttctgctagctgtgctcactgtgtgctgcagtggaacagactgtctga