

LEP-leptin Gene

Reference: PTXBC069323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEPD; OBS; leptin (murine obesity homolog); leptin (obesity homolog, mouse); obese protein; obese, mouse, homolog of; obesity factor		

For research use only.

GENE DESCRIPTION	leptin
GENE ID	3952
ORIGIN SPECIES	Human
ACCESSIONS	BC069323
VECTOR	pDONR223
NCBI SYMBOL	LEP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcattggggaacctgtgcggattcttggctttggccctatctttctatgcccaagctgtgccatccaaaagtccaagatgacacaaaacctcatcaagaattgtcaccaggatcaatgacatttcacacacgagtcagtcctccaaacagaaagtcaccgggttggaattc
attcctgggctccacccatcctgacctatccaagatggaccagacactggcagctctaccaacagatcctcaccagatgcttccagaaacgtgatccaaatccaacgacctggagaacctccgggatcttcttcacgtgctggccttcttaagagctgccacttgcctgggcccagtg
gcctggagaccttggacagcctgggggtgtcctggaagcttcaggctactccacagaggtgggtggccctgagcaggctgcagggtctctgcaggacatgctgtggcagctggacctcagccctgggtgctga