

RHOXF1-Rhox homeobox family, member 1 Gene

Reference: PTXBC069324

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OTEX; PEPP1; rhox homeobox family member 1; PEPP subfamily gene 1; ovary-, testis- and epididymis-expressed gene protein; paired-like homeobox protein PEPP-1		

For research use only.

GENE DESCRIPTION	Rhox homeobox family, member 1
GENE ID	158800
ORIGIN SPECIES	Human
ACCESSIONS	BC069324
VECTOR	pDONR223
NCBI SYMBOL	RHOXF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgttcgctcgtccagacacccgtgttctactgcctgagtgatgtatccaggtataaataagccccacacctcagctgggggcagcatcaagcgcagaaggccatgttgccaaggagctccaggcctcatgggtaatatgaacctgagggcggtgtgaaccacgagaacggcat
gaaccgcgatggcgccatgatccccgagggcggttggaaccaggagcctcggcagcagccgcagccccgcggaggagccggccaggcgcctatgagggtccgcagcccagagaacatgcagccacgaactcggcgacgaagttcacgctgttcaggtggagga
gctggaaagtgtttccgacacactcaataccctgatgtgcccaagaagggaacttgcgaaaacttaggtgtgactgaagacaagtgcgggtttgtttaagaataaaaggccagatgtaggcgacatcagagagaattaatgctcgccaatgaactacgtgtgaccagacgac
tgtgtctacatcgtcgtggactag