

KLK4-kallikrein-related peptidase 4 Gene

Reference: PTXBC069325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AI2A1; ARM1; EMSP; EMSP1; KLK-L1; PRSS17; PSTS; kallikrein; kallikrein-4; androgen-regulated message 1; enamel matrix serine protease 1; enamel matrix serine proteinase 1; kallikrein-like protein 1; prostase; serine protease 17; kallikrein related peptida		

For research use only.

GENE DESCRIPTION	kallikrein-related peptidase 4
GENE ID	9622
ORIGIN SPECIES	Human
ACCESSIONS	BC069325
VECTOR	pDONR223
NCBI SYMBOL	KLK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacagcaggaaatccctggggctgttctctgggtacatcttgggtgcgaggatcgtctctgttagctgcagccaaatcataaacggcgaggactgcagcccgactcgcagccctggcaggcggcactggtcatggaaaacgaattgtctgctcggcgctcgtg
gtgcatccgcagtggtgtgtgtcagccgcacactgttccagaactctacacatcgggctgggctgcacagcttggagccgaccaagagccagggaagccagatggtgaggccagcctctccgtacggcaccagagtaacaagaccctgtcgtcgaacacatgctcat
caagttggacgaatccgtgtccgagctgacacatccggagatcagcattgcttcgagtgccctaccgcccggggaactctgctcgtttctgctggggtctgctgagcgaacggcagaatgcctaccgtgctgagtgctgaacgtgtcgtgtgtgtgaggggtctgagtaag
ctctatgacccgctgtaccacccagcatgttctgcgccggcgagggaagaccagaaggactcctgcaacggtgactctggggggccctgactgcaacgggtactgcaggccctgtgtctttcgaaaagccccgtgtgccaagttggcgtgccaggtgtctacacacac
tgcaaatcactgagtgatagagaaaaccgtccaggccagttaa