

GGT1-gamma-glutamyltransferase 1 Gene

Reference: PTXBC069326

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD224; D22S672; D22S732; GGT; GGT 1; GTG; gamma-glutamyltranspeptidase 1; glutathione hydrolase 1; leukotriene-C4 hydrolase; testicular tissue protein Li 73; gamma-glutamyltransferase 1		

For research use only.

GENE DESCRIPTION	gamma-glutamyltransferase 1
GENE ID	2678
ORIGIN SPECIES	Human
ACCESSIONS	BC069326
VECTOR	pDONR223
NCBI SYMBOL	GGT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctctgagttctctgctgcccagctccggcccagatctgacgacacacacacccgatctctactacaagcccgagttctacacgccggatgacgggggcaactgtcacctgtctgtctgcaaggagcagtcgtgtgtccgaccagcaccatcaacctctacttgg
ctccaagggtccgctccccggtcagcgggatctgttcaataatgaaatggacgacttcagctctccagcatcaccaacgagtttgggtacccccctcacctgccaatctccagccagggaagcagccgctctctccatgtgccgacgatctgtgtggccaggacggccaggt
ccggatggtgtgtgggagctgtgtggggcacacagatcaccacggccactgcactggccatcatctacaacctctggttcggtatgacgtgaagcggaccgtggaggagccccggctgcacaaccagcttctgcccaacgtcacgacagtgagagaaacattgaccaggcagtgac
tgacgcccctggagacccggcaccatcacaccagatcgctccaccttcacgtgtgtgtgcaagccatcgctccgacggctggtgtgtggcagctgcctcggactccaggaaaggcggggagcctgcccgtactga