

IFNA21-interferon, alpha 21 Gene

Reference: PTXBC069329

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFN-alpha1; LeIF F; leIF-F; interferon alpha-21; IFN-alpha-21; interferon alpha-F; leukocyte interferon protein; interferon alpha 21		

For research use only.

GENE DESCRIPTION	interferon, alpha 21
GENE ID	3452
ORIGIN SPECIES	Human
ACCESSIONS	BC069329
VECTOR	pDONR223
NCBI SYMBOL	IFNA21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgtcctttcttactgatggccgtgctggtgctcagctacaaatccatctgttctctgggctgtgatctgcctcagaccacagcctgggtaataggaggccttgatactcctggcacaatgggaagaatctctcttctcctgcctgaaggacagacatgactttggattcccca
ggaggagtgtgatggcaaccagttccagaaggctcaagccatctctgtctccatgagatgatccagcagaccttcaatctcttcagcacaaggactcatctgctacttgggaacagagcctcctagaaaaatttccactgaactaaccagcagctgaatgacctgggaagcctgctgata
caggaggttggggtggaagagactccccgatgaatgtggactccatcctggctgtgaagaatacttccaaagaatcactcttctgacagagaagaatacagccctgtgctctgggaggtgtgacagcagaatcatgagatcctctctttatcaaaaattttcaagaagattaagg
aggaaggaatga