

## FGF23-fibroblast growth factor 23 Gene

Reference: PTXBC069333

|                |                                                                                                                           |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ADHR; FGFN; HPDR2; HYPF; PHPTC; fibroblast growth factor 23; phosphatonin; tumor-derived hypophosphatemia inducing factor |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | fibroblast growth factor 23                                                                       |
| <b>GENE ID</b>          | 8074                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC069333                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FGF23                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgttggggcccgctcaggtctgtgctgtgctgagcgtctgcagcatgagcgtctcagagcctccaatgcctcccaactgctcggtccagctgggtggcctgatccactgtacacagccacagcaggaacagctaccacctgcagatccacaagaatggccat  
gtggatggcgcaccccatcagaccatctacagtgcctgatgatcagatcagaggtgctgctgtgtgattacaggtgtgatgagcagaagatacctctgcatggatttcagaggcaacattttggatcacactatttgcacccggagaaactgcaggttccaacaccagacgtggaa  
aacgggtacgacgtctaccactctcctcagatcacttctgtgctgctggcggcggaagagagccttctgcccagcatgaaccacccccgactccagttctgtccggagggaacgagatcccccttaattcactcaacacccccataccacggcgacacccggagcgcc  
gaggacgactcggagcgggacccccgaacgtgctgaagccccggcggcgatgacccccggcctctgttcacagagctcccgagcgccgaggacaacagcccgatggccagtgacccattaggggtgtcagggcggtcagtgaaacacgacgtggggaa  
cgggccccggaaggctgcgcccccttgcgaagtcatctag