



Reference: PTXBC069345

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAS-L; MRG; dJ994E9.2; mas-related G-protein coupled receptor MRG; MAS-R; MAS1 oncogene-like; MAS1-like; mas-related G protein-coupled MRG; MAS1 proto-oncogene like, G protein-coupled receptor		

For research use only.

GENE DESCRIPTION	MAS1 oncogene-like
GENE ID	116511
ORIGIN SPECIES	Human
ACCESSIONS	BC069345
VECTOR	pDONR223
NCBI SYMBOL	MAS1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]