



Reference: PTXBC069351

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARCI5; INLNE; LI3; cytochrome P450 4F22; cytochrome P450, family 2, subfamily E, polypeptide 2 homolog; cytochrome P450, family 4, subfamily F, polypeptide 22; cytochrome P450 family 4 subfamily F member 22		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 4, subfamily F, polypeptide 22
GENE ID	126410
ORIGIN SPECIES	Human
ACCESSIONS	BC069351
VECTOR	pDONR223
NCBI SYMBOL	CYP4F22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]