



Reference: PTXBC069356

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	5-HT2A; HTR2; 5-hydroxytryptamine receptor 2A; 5-HT2 receptor; 5-hydroxytryptamine (serotonin) receptor 2A, G protein-coupled; serotonin 5-HT-2A receptor		

For research use only.

GENE DESCRIPTION	5-hydroxytryptamine (serotonin) receptor 2A
GENE ID	3356
ORIGIN SPECIES	Human
ACCESSIONS	BC069356
VECTOR	pDONR223
NCBI SYMBOL	HTR2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr