

XCL2-chemokine (C motif) ligand 2 Gene

Reference: PTXBC069360

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCM-1b; SCM1B; SCYC2; cytokine SCM-1 beta; XC chemokine ligand 2; c motif chemokine 2; chemokine (C motif) ligand 2; small inducible cytokine subfamily C, member 2; X-C motif chemokine ligand 2		

For research use only.

GENE DESCRIPTION	chemokine (C motif) ligand 2
GENE ID	6846
ORIGIN SPECIES	Human
ACCESSIONS	BC069360
VECTOR	pDONR223
NCBI SYMBOL	XCL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagactctcatctggccctccttggcatctgctctcactgcatactgtggaaggtgtgaggagtggaagtcacataggagacgtgtgagcctcactaccagcgactgccagtttagcagaatcaagacctacaccatcacggaaggctccttgagagcagtaattttattac
caaacgtggcctaaaagtctgtgctgatccacaagccacgtgggtgagagacgtggtcaggagcatggacaggaaatccaacaccagaaataacatgatccagaccaagccaaggaaccagcaatcgaccaatacagctgtgaccctgactggctag