

GNRH2-gonadotropin-releasing hormone 2 Gene

Reference: PTXBC069362

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GnRH-II; LH-RHII; progonadoliberin-2; GnRH-associated peptide 2; GnRH-associated peptide II; Gonadoliberin II; Gonadoliberin-2; luliberin II; luteinizing hormone-releasing hormone II; progonadoliberin II; gonadotropin releasing hormone 2		

For research use only.

GENE DESCRIPTION	gonadotropin-releasing hormone 2
GENE ID	2797
ORIGIN SPECIES	Human
ACCESSIONS	BC069362
VECTOR	pDONR223
NCBI SYMBOL	GNRH2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagctccagcgaggcctcctgctcctgctgctgctgactgccaccttgaccctcagaggetcagcactggctccatggctgtaccctggaggaagcagagccctcagctcagccagatccccaatgcccttaggccccaggaaggccctggacactgcagca
ggcagccccagtcagactgccatggcctcccaagtgtgcctggctccctggacgacagcatgccctggaggcaggaccacggccagtggtcccttcacaggaagcagacactggcagggacactgctgaccgcagcccgagagccccgcccgccccgcatcctcca
ataaagtgtga