

CAV3-caveolin 3 Gene

Reference: PTXBC069368

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LGMD1C; LQT9; VIP-21; VIP21; caveolin-3; M-caveolin; cavolin 3		

For research use only.

GENE DESCRIPTION	caveolin 3
GENE ID	859
ORIGIN SPECIES	Human
ACCESSIONS	BC069368
VECTOR	pDONR223
NCBI SYMBOL	CAV3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggcagaagacacagatctcagggccagatcgtcaaggatccactgcaaggagattgacctggtgaaccgagaccccaagaacattaatgaggacatagtcagggtgattttgaagacgtgatcgcagagcctgtgggcacctacagcttgacggcgtgtggaaggt
gagctacaccaccttcactgtctccaagtactggtgtaccgtctgtgtccacgtctgtggcgctccactggccctgctctggggcttctgttcgctgcactctcttctgccacatctgggcgggtggtgccatgcaataagagctacctgatcgagatccagtcgcatcagccacatctactc
actctgcacccgacactcttgcacccactcttcgcgccctggggccaggtctgcagcagcatcaagggtggtgctgcggaaggagggtctaa