

TAS2R13-taste receptor, type 2, member 13 Gene

Reference: PTXBC069376

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T2R13; TRB3; taste receptor type 2 member 13; taste receptor, family B, member 3; taste receptor, type 2, member 13; taste 2 receptor member 13		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 13
GENE ID	50838
ORIGIN SPECIES	Human
ACCESSIONS	BC069376
VECTOR	pDONR223
NCBI SYMBOL	TAS2R13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaagtgcctgcccagatcttctcactcttgaataattgcagaattcataattgggaatttgagcaatggattatagtagtcatcaactgcattgactgggtcagtaaaagagagctgtctcagtcgataaactcctcattatctggcaatccagaattgggctgatctgggaaatattag
taagtgggttttagctctgcattatctagccatatttggctggaacaggattaagaattatgatttttagctggatagtttctaatcacttcaatctctggctgtacatacttcagcatctttatttgcataaatagcgagtttcttagccctgcttttctatttgaagtggagagtaaacaaagtga
ttctgatgatactgctaggaaccttggcttcttattttaaactctgatacaataaacaatgcataaaaagactggctggaccgatatgaagaacacacaaacttgaatttcagtagtgacttgaacattttcagtgctgggtcaaaatcactatgactatgttcagtcctaacaccatttacttgggc
cttcacatcttttctctgttaattttctccctgcagaacatctccagaaaatgcaactcaattacaaaggacacagagacccaggaccaaggtccatacaaatgcttgaattgtgatctcattctttattctatgctagttttctatgtgttctcatatcatggatttctgagctgtatcagaa
cacagtgatctacatgctttgtgagacgattggagcttctctcttcaagccactcttctctgattctaggaacgctaagtaagacagggcctttcttttgggtggcagctaaggtatgggctaaacgatga