

SPANXA1-sperm protein associated with the nucleus, X-linked, family member A1 Gene

Reference: PTXBC069393

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT11.1; NAP-X; SPAN-X; SPAN-Xa; SPAN-Xb; SPANX; SPANX-A; sperm protein associated with the nucleus on the X chromosome A; SPANX family member A; SPANX family, member A1; cancer/testis antigen 11.1; cancer/testis antigen family 11, member 1; nuclear-associ		

For research use only.

GENE DESCRIPTION	sperm protein associated with the nucleus, X-linked, family member A1
GENE ID	30014
ORIGIN SPECIES	Human
ACCESSIONS	BC069393
VECTOR	pDONR223
NCBI SYMBOL	SPANXA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaacaatccagtgcggcggtgaagaggagcgtccctgtgattccaacgagccacgagatgatccggagaccccaactggggactcagaccgcaacctgtcctaaaaaatgaaacatctgagtcctcgaccatactagtgttcgctacaggaggaactttaa
aagaacatctccagaggaactgttgaatgaccacgccgagagaacagaatcaacccctccaatggaggagggaattcatggaataatggttgaatacctgcaaagtag