



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

P53AIP1-p53-regulated apoptosis-inducing protein 1 Gene

Reference: PTXBC069399

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P53AIP1; p53-regulated apoptosis-inducing protein 1; tumor protein p53 regulated apoptosis inducing protein 1		

For research use only.

GENE DESCRIPTION	p53-regulated apoptosis-inducing protein 1
GENE ID	63970
ORIGIN SPECIES	Human
ACCESSIONS	BC069399
VECTOR	pDONR223
NCBI SYMBOL	P53AIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagatctctctcgtgaggtgcagctcagatctgctcagagctctctcagtggtggggcaggaggcaggggcctgggcaggggagaccagaacctctcgtgtagctccctggaatggcaggggctcagacacacacacctgctggglaagccctgcagtgaaaccgagacggctttggcc
tgcacagaccccggcagagctctgctctaccgtgtgtccgacatcccaagttttcagatcaccaggaccagctgcacagcatctgggtcctcagacgtcgatcgcgactctcccaagtcagacattggacagagctctgggaactga