

CRYBA4-crystallin, beta A4 Gene

Reference: PTXBC069404

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTRCT23; CYRBA4; MCOPCT4; beta-crystallin A4; beta crystallin A4 chain transcript PS; beta crystallin alpha 4 chain; beta-A4 crystallin; crystallin, beta polypeptide A4; eye lens structural protein; crystallin beta A4		

For research use only.

GENE DESCRIPTION	crystallin, beta A4
GENE ID	1413
ORIGIN SPECIES	Human
ACCESSIONS	BC069404
VECTOR	pDONR223
NCBI SYMBOL	CRYBA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgcaatgcacaaagtcagcgggacccctggaagatgggtgtgtgggatgaggacggctccaggccggcgccaggtcacggccgagtgccccagcgtgctggagcttgcttcgagactgtgcgatcttgaagtgtgagtgagcgtgggtggcttcgagcatg
ctggcttcaaggcgacgagtaacttcggaacgaggcgaatatccaagctgggatgctggggcggaacacggcctaccccgccgagaggtcacctccttcggcctgcgcctgtgctaaccacgtgactcgaggctgacaatcttcgagcaagagaactctgggcaagaa
aggagagctgagcgtatgactatcctccctccaggccatgggatgggaaggcaatgaagtagggctcctccagctccactctggggcctgggtttgctccagtttcgggctaccgaggatttcagtatgtgctggaatgcgacacattccgggtgactacaacattccgggagtggg
gctctcatgccccgacctccaggtgcagagcatccgaggtatccagtcagtga