

FAM12A-family with sequence similarity 12, member A Gene

Reference: PTXBC069407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM12A; EP3A; HE3-ALPHA; HE3A; HE3ALPHA; RAM1; epididymal secretory protein E3-alpha; epididymis-specific 3 alpha; family with sequence similarity 12, member A; human epididymis-specific 3 alpha; human epididymis-specific protein 3-alpha; ribonuclease A M		

For research use only.

GENE DESCRIPTION	family with sequence similarity 12, member A
GENE ID	10876
ORIGIN SPECIES	Human
ACCESSIONS	BC069407
VECTOR	pDONR223
NCBI SYMBOL	FAM12A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacatcctctctaaagattggggcactactcttggccctgtttgcatcctttgcaggctgtgtgtatagtaacaacatttactggagagaattcataaaacttcattacttaagtccaagtcgagaattcaaaagagtacaaatgtgatgtcctcatgagagaaaaaggcctctgaaaggcaa
gagctttcatatgttcatctatagcttatgggttcaaaatcagcgtgcatgcatcaatgagaaggggagtgaccgatatagaaatgcataatgtatgggccccaggtgccctcaaaagtactcgagtgctactgggagaagtacaacaatagggtacacagagagcagaagcttcagctacattgaa
ttccattgtggcgtagatggatgttgataacatagaagacctgaggattatagaacctatcagcaactag