



Reference: PTXBC069420

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIEE12; PI-PLC; PLC-154; PLC-I; PLC-beta-1; PLC154; PLCB1A; PLCB1B; 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1; 1-phosphatidyl-D-myo-inositol-4,5-bisphosphate; inositoltrisphosphohydrolase; monophosphatidylinositol phosphodiesterase;		

For research use only.

GENE DESCRIPTION	phospholipase C, beta 1 (phosphoinositide-specific)
GENE ID	23236
ORIGIN SPECIES	Human
ACCESSIONS	BC069420
VECTOR	pDONR223
NCBI SYMBOL	PLCB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]