

MT1B-metallothionein 1B Gene

Reference: PTXBC069421

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT-1B; MT-IB; MT1; MT1Q; MTP; metallothionein-1B; metallothionein 1B (functional); metallothionein 1Q; metallothionein-IB; metallothionein 1B		

For research use only.

GENE DESCRIPTION	metallothionein 1B
GENE ID	4490
ORIGIN SPECIES	Human
ACCESSIONS	BC069421
VECTOR	pDONR223
NCBI SYMBOL	MT1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcccaactgctcctgcaccacaggtggctcctgtgcctgcgccggctcctgcaagtgcgaagagtgcaaatgtacctcctgcaagaagtgcctgctcttgcctgccccgtgggctgtgccaaagtgtgccagggctgtgtctgcaaggctcatcagagaagtgccgctgctgtg
cctga