

DEFA1-defensin, alpha 1 Gene

Reference: PTXBC069423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEF1; DEFA2; HNP-1; HP-1; HP1; MRS; neutrophil defensin 1; defensin, alpha 1, myeloid-related sequence; defensin, alpha 2; myeloid-related sequence; defensin alpha 1		

For research use only.

GENE DESCRIPTION	defensin, alpha 1
GENE ID	1667
ORIGIN SPECIES	Human
ACCESSIONS	BC069423
VECTOR	pDONR223
NCBI SYMBOL	DEFA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaccctcgccatcctgctgccattctctgtggccctgcaggcccagctgagccactccaggcaagagctgatgaggttgcagccccggagcagattgcagcggacatcccagaagtgtttccctgcatggacgaaagcttgctccaaagcatccaggctcaaggaaaaacatggcctgctattgcagaataccagcgtgcattgcaggagaacgtcctatggaacctgcatctaccagggaaagactctggcattctgctctga