

IFNA13-interferon, alpha 13 Gene

Reference: PTXBC069427

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	interferon alpha-1/13; IFN-alpha-1/13; interferon alpha-D; leIF D; interferon alpha 13		
For research use only.			

GENE DESCRIPTION	interferon, alpha 13
GENE ID	3447
ORIGIN SPECIES	Human
ACCESSIONS	BC069427
VECTOR	pDONR223
NCBI SYMBOL	IFNA13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgcgctttgctttactgatggcctgggtgctcagctgcaagtcagctctctgggctgtgatctccctgagaccacagcttgataacaggaggaccttgatctcctggcacaaatgagcagaatctctctctctctgtctgatggacagacatgacttggattcccc
aggaggagtttgatggcaaccagttccagaaggctccagccatctctgtctccatgatgctgacccagcagatcttcaacctcttaccacaaaagattcatctgctctgggatgaggacctctagacaaattctgcaccgaactctaccagcagctgaatgacttggagcctgtgtgatg
caggaggagagggtgggagaaactccctgatgaatgcggactccatcttggctgtgagaataacttccgaagaatcactctctatctgacagagaagaatacagcccttgctcctgggaggtgtcagagcagaatcatgagatccctctttatcaacaacttgcagaagatta
aggagggaaggaataa