

KRTAP3-3-keratin associated protein 3-3 Gene

Reference: PTXBC069448

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP3.3; KRTAP3.3; keratin-associated protein 3-3; high sulfur keratin-associated protein 3.3; keratin-associated protein 3.3; keratin associated protein 3-3		

For research use only.

GENE DESCRIPTION	keratin associated protein 3-3
GENE ID	85293
ORIGIN SPECIES	Human
ACCESSIONS	BC069448
VECTOR	pDONR223
NCBI SYMBOL	KRTAP3-3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattgctgtcctctcagagctgcagtgctccaccggcctgccaccacatctgctctccgacaaatctgtgctgtggagctgtcctgccagcactgccacacacagtttggttaactggagccacctgtgtgacaaactgtccccaccctgccacattctcagccctgtgtgcccaactgtcttctgtctcaactctctccagccaactccagcctggagacctcaacctcaccacttcaactcagccctgtgtgagccctgctctccaagagctgctaa