

RP11-265F14.2-elastase 2B Gene

Reference: PTXBC069455

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RP11-265F14.2; ELA2B; chymotrypsin-like elastase family member 2B; elastase 2B; pancreatic elastase IIB (ELA2B); chymotrypsin like elastase family member 2B		

For research use only.

GENE DESCRIPTION	elastase 2B
GENE ID	51032
ORIGIN SPECIES	Human
ACCESSIONS	BC069455
VECTOR	pDONR223
NCBI SYMBOL	RP11-265F14.2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattaggaccctgctgctccacttgggtggagccctcagttgtgggtctccacttacgcgcctgatatgtctagatgctggaggtgaagaagcagagcccaacagctggccctggcaggtctccctgcagtacagctccaatggccagtggtaccacacctgcggagggt
ccctgatagcccaacagctgggtctgacggctgcccactgcatcagctcctccgggatctaccgcgtgatgctgggcccagcacaacctctacgttcagagtcgggctcgtggccgtcagtgctcttaagattggtgacaaagactggaactccgaccaggtctccaaagggaacg
acattgcccgtctcaactggctaaacccgtctccctcaccgacaagatccagctggcctgctcctcctgcccggcaccattctacccaacaactaccctgctacgtcacgggtggggaaggctgcagaccaacgggctctccctgatgacctgaagcaggccggttctgtgtgt
ggactatgccacgtgctcagctctggctgtggggcagcaccgtgaagacgaatatgatctgtgctgggggtgatggcgtgatatgcacctgcaacggagactccgggtggccgtgaactgtcaggcatctgacggccgtgggaggtgcatggcatggcagcctcacgtcggtc
cttggttgaactactactacaagccctccatttcacgggtctccaactacaacgactggatcaattcgggtgattgcaataactaa