

GIMAP6-GTPase, IMAP family member 6 Gene

Reference: PTXBC069461

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IAN-2; IAN-6; IAN2; IAN6; GTPase IMAP family member 6; immune associated nucleotide 2; immune associated nucleotide 6; immune-associated nucleotide-binding protein 6; GTPase, IMAP family member 6		

For research use only.

GENE DESCRIPTION	GTPase, IMAP family member 6
GENE ID	474344
ORIGIN SPECIES	Human
ACCESSIONS	BC069461
VECTOR	pDONR223
NCBI SYMBOL	GIMAP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaagaagaatgaacaaattcccaggagaatccccagaagagctgtccagatcctgtgctggagctgtcaggaggtctaaggagaaagaacagaagacccaaggagactgaggtcattctcatgggaaaacaggagtggaagagtgcaacaggaaaca
gcatctcggcaggagactctcgaagtctaaacacagcaccagaccgtgacaaagacctccagagacggagccgagagtggtgggaaaggagcttgaggtgattgacacccaacattctgtccccagggtctcgcagaggtggcagacgtatctgccaagccatcgttta
tccgccccaggccccacgcccgtgctcctggtgacacaactgggcccgttcacggatgaggaatcagcaggtggtcaggcgcctgcaggaggtcttggagtggtgggttctgggtcacaccatcctggtgttccccggaaggaaagacctggtggtggtcctcctggaagactatgtgc
gagagaccaacaacaggcccttgctgctggtgatgtgaccttgcacggcgccattgcccgttcaacaacagggcacagggggaggagcaggagccaactgcgagagctcatggagaaagttgaagccattatgtgggaaaacgaaggagattattacagcaacaaggcttacc
aatataccagcaaaactttcggctga