



Reference: PTXBC069465

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ECNOS; nitric oxide synthase, endothelial; EC-NOS; NOS type III; NOSIII; cNOS; constitutive NOS; endothelial NOS; nitric oxide synthase 3 (endothelial cell); nitric oxide synthase 3 transcript variant eNOS-delta20; nitric oxide synthase 3 transcript varia		

GENE DESCRIPTION	nitric oxide synthase 3 (endothelial cell)
GENE ID	4846
ORIGIN SPECIES	Human
ACCESSIONS	BC069465
VECTOR	pDONR223
NCBI SYMBOL	NOS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]