

HIST1H4B-histone cluster 1, H4b Gene

Reference: PTXBC069467

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4/I; H4FI; histone H4; H4 histone family, member I; histone 1, H4b; histone cluster 1, H4b; histone cluster 1 H4 family member b		

For research use only.

GENE DESCRIPTION	histone cluster 1, H4b
GENE ID	8366
ORIGIN SPECIES	Human
ACCESSIONS	BC069467
VECTOR	pDONR223
NCBI SYMBOL	HIST1H4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtgtcgcgcaaggcggtaaggtttgggtaaggaggtgccaaagcgtcaccgaaaagtgcgggataacatcaaggcatcaccaaacggccattcggcgcttgctaggcgtgtggggttaagcgaattccggttgattatgaggagactcgtggcgttctcaag
gtgtttctggagaacgtgatccgggacgcccgtgacctacacggagcagcgaagcgcaagactgtcactgcatggatgtgtttacgcgtcaagcgtcaaggagcgcactctgtacggcttcggcggttaa