

OCM-oncomodulin Gene

Reference: PTXBC069468

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OCM1; ONCM; oncomodulin-1; beta parvalbumin; hCG18255; parvalbumin beta; oncomodulin		

For research use only.

GENE DESCRIPTION	oncomodulin
GENE ID	654231
ORIGIN SPECIES	Human
ACCESSIONS	BC069468
VECTOR	pDONR223
NCBI SYMBOL	OCM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcatcacgagtgctcagtgctgacgacattgcagcagcgtccaggaatgccgagaccagacactttgaaccccaaaattctccagacatcaggcctctccaagatgtcagccaatcagggtgaaggatgtttccggttcatagacaacgaccagagcgggtacctggatg
aagaagagcttaagtgtttctccagaagtgtgagagtggtgccagagaactgaccgagtcagaaaccaagtccttgatggctgcggcgataatgatggagatgggaaattggagcagaggaaatccaggaaatgggtgcattcttaa