

HIST1H1A-histone cluster 1, H1a Gene

Reference: PTXBC069492

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1.1; H1A; H1F1; HIST1; histone H1.1; H1 histone family, member 1; histone 1, H1a; histone H1a; histone cluster 1, H1a; histone cluster 1 H1 family member a		

For research use only.

GENE DESCRIPTION	histone cluster 1, H1a
GENE ID	3024
ORIGIN SPECIES	Human
ACCESSIONS	BC069492
VECTOR	pDONR223
NCBI SYMBOL	HIST1H1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgaaacagtgctccgcccccgcctctctgctgctgagaaacctttagctggcaagaaggcaagaacctgtaaggctgcagcagctccaagaaaaacccgtggcccttcgtgcagagctgatcgtgcaggctgtcttcctcctaaggagcgtggtgtgtgt
cgftggcagctcttaaaaaggcgtggcgccgcaggctacgacgtggagagaacaacagccgcattaaagctgggcattaaagagcctgtaagcaagggaacgttggtgcagacaagggtaccggagcctcggttcctcaagctcaacaagaaggcgtcctccgtggaacca
agcccgccgctcaaaagtgctacaaaaactaaggcaacgggtgcatctaaaaagctcaaaaaggccacggggctagcaaaaagagcgtcaagactccgaaaaaggctaaaaagcctgaggcaacaaggaaatcctcaagaatccaaaaaacccaaactgtaaaagccca
gaaagtagctaaaaagccctgctaaagctaaagctgtaaaacccaaggcgcccaaggctagggtgacgaagccaagactgccaaaccaagaaagcggcacccaagaaaaagtaa