

PRSS3-protease, serine, 3 Gene

Reference: PTXBC069494

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MTG; PRSS4; TRY3; TRY4; trypsin-3; brain trypsinogen; mesotrypsin; mesotrypsinogen; pancreatic trypsinogen III; protease, serine, 4 (trypsin 4, brain); trypsin III; trypsin IV; trypsinogen 4; trypsinogen 5; trypsinogen IV; protease, serine 3		

For research use only.

GENE DESCRIPTION	protease, serine, 3
GENE ID	5646
ORIGIN SPECIES	Human
ACCESSIONS	BC069494
VECTOR	pDONR223
NCBI SYMBOL	PRSS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatccattctgatccttgcctttgtggagctgctgttctgtcccccttgacgatgacaagattgtgggggtacacctgtgaggagaattctctccccaccaggtgtccctgaattctggctcccacttctgcgggtggtccctcatcagcgaacagtggtggtatcagcagctc
actgctacaagaccgcatccagggtgagactgggagagcacaacatcaaaagtcctggagggaatgagcagttcatcaatgcgccaagatcatccgccaccctaaatacaacaggagacactctggacaatgacatcatgctgatcaaaactctctcacctgccgtcatcaatgcccgcg
tgtccaccatctctgtcccaccgccccctccagctgctggcactgagtgccctcatctccgctggggcacaactctgagcttgggtgctgactaccagacgagctgaagtgcctggatgctccggtgctgaccagagctgagtgtaaaagcctcctaccctggaaaattaccaacagcatg
ttctgtgtgggtctcttgaggagggaaggattctgccagcgtgactctgggtggcctgtggtctgcaacggacagctccaaggagttgtctctggggccatggctgtgctggaagaacaggcctggagtctacccaaggtctacaactatgtggactggattaaggacaccatcgc
tgccaacagctaa