

CST8-cystatin 8 (cystatin-related epididymal specific) Gene

Reference: PTXBC069496

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRES; CTES5; cystatin-8; cystatin 8 (cystatin-related epididymal specific); cystatin-related epididymal spermatogenic protein; cystatin 8		

For research use only.

GENE DESCRIPTION	cystatin 8 (cystatin-related epididymal specific)
GENE ID	10047
ORIGIN SPECIES	Human
ACCESSIONS	BC069496
VECTOR	pDONR223
NCBI SYMBOL	CST8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgccaggtgccgtgctctccctgatcctcaccattccctggccctgggccaggaagacccaaaaagaatgagcgggggtgctgaggaattaaaccctcaatgcctcaaatccaacgtgaagcagtgctgtgtttgccatgcaagaatacaaaaagagc
gaggacaagtatgtctcctgggtgcaagacactgcaagcccagctcaggtcacaatctctggaatacctattgatgtagaattgcccgcagcgattgcagaagcctttaagcactaatgaaatctcgccatcaagaaaactccaagctgaaaaggaaataagctgcagctttt
gtaggagcactccctggaatgggaattcactgtgatggagaaaaagtgaagatgcttaa
```