

HOXA6-homeobox A6 Gene

Reference: PTXBC069497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein HOXA6; HOX1; HOX1.2; HOX1B; homeobox protein Hox-A6; homeo box 1B; homeo box A6; homeobox protein Hox-1B; homeobox A6		

For research use only.

GENE DESCRIPTION	homeobox A6
GENE ID	3203
ORIGIN SPECIES	Human
ACCESSIONS	BC069497
VECTOR	pDONR223
NCBI SYMBOL	HOXA6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttcctattttgtaaccactttccccgggagcctccagcggcaggactcctcttggccagctgccccctaccagctggctatgacgcgtgagcccttccccggcctcgtacggggcgtcgtcgtcctccggacaagacgtacacctcacttgttttaccacagtcac
actcggctcgtgctgcaaccggggtcctacgagtagggcctcgtgtttctattctgataaggacacagtgccgctcgcctcggcagtggaagcagagggcccccgggactacctgacatttttccccgagcagcagtaaaacccgacagcagcggcgaggca
aagcactccatgacgaaggcggccgaccggaagtacacgagcccggtttacccttggaatgcagcggatgaactcctgcgcgggtgctgtgtatgggagccatgggcggccgagggccgagacacacgcgtaccagacactggagctggagaaggagttcacttcaaccgctac
ctgacacggcggccgcatcgagatgccaacgcgtctgctcaccgagcgcagatcaagatctggttcagaaccgcccgaatgaagtggaaaaaggaaaaacaagctcatattccacgcagcccgaggaggactcagaggcaaggcggcgagtag