

MDS1-myelodysplasia syndrome 1 Gene

Reference: PTXBC069498

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDS1 and EVI1 complex locus; MDS1 and EVI1 complex locus protein MDS1; MDS1 and EVI1 complex locus protein EVI1; MDS1-EVI1; MDS1; AML1-EVI-1; EVI1; KMT8E; PRDM3; RUSAT2; AML1-EVI-1 fusion protein; PR domain 3; ecotropic virus integration site 1 protein ho		

For research use only.

GENE DESCRIPTION	myelodysplasia syndrome 1
GENE ID	4197
ORIGIN SPECIES	Human
ACCESSIONS	BC069498
VECTOR	pDONR223
NCBI SYMBOL	MDS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagatccaaaggcagggaaggaaactggccacaataatgagtggtatggaactaccctgaaataaccttggagaagatccagatgcagatggagtagccagcactccctccctcaatattcaagagccatgctctcctgccatccagtgaagcattcactccaaaggagggttctccttacaagcccccatctacatccctgatgataccccatctcctgagtttgaactcgagagtcacaatagcctggggcaggactaggaatatggaccaaagggaagatcgaagtaggtgaaaagtgtggcccttatgtgggagagcagaggtcaaacctgaaagacccagttatggatgggaggtacatctccaaggtctcggagggtgaagcgttcactctgtgtgtatgtgggaagagaagctcagacgtaggaatagcctctcagagctgatgtctacatgcctggactgcagtggtgccttcctctcgtag