

BMF-Bcl2 modifying factor Gene

Reference: PTXBC069505

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bcl-2-modifying factor; Bcl2 modifying factor		
For research use only.			

GENE DESCRIPTION	Bcl2 modifying factor
GENE ID	90427
ORIGIN SPECIES	Human
ACCESSIONS	BC069505
VECTOR	pDONR223
NCBI SYMBOL	BMF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagccatctcagtggtggaggagctggaggatgatgtttccaaccagaggatggggagccggtagcccaacccgggagcttgcctctgtgacctgtttgccagagcctactggagctccccctcagccgacttcagctctccctcaccactgctgtggccctggccttcg
accacccagccaggaagacaagctaccagactctcagccagcctccccagccaggtgtcatgtgccttggggtagctgaggaaacccagcagctttttatggcaatgctggctatcggttctctccctgccagtttccagcagcttgccttggggagcagccccc
gaaggcagtgaggcaatcaagcagaggtacagattgcccgaagctcagtcattgcagaccagttccacggcttcatgtgcagcaacaccagcagaacaaaatcgtgtgtgggcagatcctcctcttcctgcacaacctgtcttgaatggagaagagaacaggaggggc
aggccctaggtga