

CST5-cystatin D Gene

Reference: PTXBC069514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cystatin-D; cystatin 5; cysteine-proteinase inhibitor; cystatin D		

For research use only.

GENE DESCRIPTION	cystatin D
GENE ID	1473
ORIGIN SPECIES	Human
ACCESSIONS	BC069514
VECTOR	pDONR223
NCBI SYMBOL	CST5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtggcccatgcacccccactgctgctgactgcttgatggtggccgtggccgggagtgccctggcccaatctaggacctggcaggtggcatccatgccacagacctcaatgacaagagtgtgcagtggtccctggactttgccatcagcaggtacaacaagggtcattaata
aggatgagtactacagccgacctctgcaggtgatggtgctgctaccagcagatcgtgggtgggtgaactactacttcaatgtgaagttcggtcgaaccacatgcaccaagtcaccagcccaactggacaactgtcccttcaatgaccagccaaaactgaaagaggaagagttctctcttc
cagatcaatgaagttccctgggaggataaaatttcattctgaactacaagtcggaagtag