

HESX1-HESX homeobox 1 Gene

Reference: PTXBC069515

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANF; CPHD5; RPX; homeobox expressed in ES cells 1; Rathke pouch homeobox; hAnf; homeobox protein ANF; HESX homeobox 1		

For research use only.

GENE DESCRIPTION	HESX homeobox 1
GENE ID	8820
ORIGIN SPECIES	Human
ACCESSIONS	BC069515
VECTOR	pDONR223
NCBI SYMBOL	HESX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctccagcctcaggagcgtcagctcggggaacaaaccccaactgtctctttcaattgagagaatcttaggactggaccagaagaagactgtgttcattatgaaacccacaggccctgggcagacacctgcagctcatcagggaagatggttaactatgtctacat
gtcccaaatctccagtggtgattcattccctagcgtgtggtgatcaccaatgccagaagaagagcttcgaaatgaaaaattactttcagcctcagaagactgtctttgaaaagagagttgagttggtatagagccgaagaccaaagactcttttactcaaaacagattgaagtgtta
gaaaatgtctttagagtaaacgtatcctggtatcgaatattagagaagacttagctcaaaaattgaatctagaggagacagaatccagattggtttcaaaatcggcgtgcaaaactgaaaaggtcccatagagaatcacagtttctaattggcgaaaaaaattcaacacaaatctgctggaa
tag