

HIST1H1T-histone cluster 1, H1t Gene

Reference: PTXBC069517

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1.6; H1FT; dJ221C16.2; histone H1t; H1 histone family, member T (testis-specific); histone 1, H1t; histone cluster 1, H1t; testicular H1 histone; histone cluster 1 H1 family member t		

For research use only.

GENE DESCRIPTION	histone cluster 1, H1t
GENE ID	3010
ORIGIN SPECIES	Human
ACCESSIONS	BC069517
VECTOR	pDONR223
NCBI SYMBOL	HIST1H1T
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttttcgctggttattacatcttgcgtttctctgtttatgtctgaaccgtgcctgcagcttgcagctgctggtgtagccgctatggagaaactccaaccaagaagcaggagggaagccgctggtgataagtgcgaagtgcgaagtgcgaacctctgtgtccaagttgatca
ccgaggcccttccagtgtcacaggaaacagtaggtatgtcttgggttcgctcaagaagcattggcctgctggtgctacgacgtagagaagaataacagccgcatcaaaactgctcctcaagagcttagtgaacaagggaatcctggtgcaaacagggtgactgtgtctccggttccttt
aagcttagtaagaaggtgattcctaatactaccagaagcaggctaaaaagtcagtttctgccaagaccaagaagctgttttatccagggactccaagtcacaaagactgctaaaaccaataagagagccaagaagccgagagcgacaactcctaaaactgttaggagcgggagaaag
gctaaaggagccaagggttaagcaacagcagaagagccagtgaaaggcaggcttcgaagtcaaaattgacccaacatcatgaagttaatgttagaaagccacatctaagaagtaa