

CCL28-chemokine (C-C motif) ligand 28 Gene

Reference: PTXBC069532

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CC chemokine CCL28; CCK1; MEC; SCYA28; C-C motif chemokine 28; chemokine (C-C motif) ligand 28 splice variant chi; mucosae-associated epithelial chemokine; small inducible cytokine A28; small inducible cytokine subfamily A (Cys-Cys), member 28; C-C motif		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 28
GENE ID	56477
ORIGIN SPECIES	Human
ACCESSIONS	BC069532
VECTOR	pDONR223
NCBI SYMBOL	CCL28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcagagaggactcgccatcgctggccttgctgtgtgcggccctacatgctcagaagccatactccattgctccagctgttgacggaggtttcacatcatattccagaaggtcctggaaagagtgaatatgtgtgcatccagagagctgatggggattgtgacttggtgctgctgtcatccttcatgtcaagcgcagaagaatctgtgtcagcccgcaaccatactgttaagcagtggaatgaaagtgcaagctgccaaagaaaatgtgtaaggaaatgtttgccacaggagaacacacattggcaagaggaaacagtaacagggcacatcaggggaaacacgaacatacggccataaaactccttattag