



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

UPK3A-uroplakin 3A Gene

Reference: PTXBC069544

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

UP3A; UPIII; UPIIIA; UPK3; uroplakin-3a; uroplakin 3; uroplakin III; uroplakin 3A

For research use only.

GENE DESCRIPTION	uroplakin 3A
GENE ID	7380
ORIGIN SPECIES	Human
ACCESSIONS	BC069544
VECTOR	pDONR223
NCBI SYMBOL	UPK3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1