



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IL25-interleukin 25 Gene

Reference: PTXBC069565

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

IL17E; interleukin-25; interleukin-17E; interleukin 25

For research use only.

GENE DESCRIPTION	interleukin 25
GENE ID	64806
ORIGIN SPECIES	Human
ACCESSIONS	BC069565
VECTOR	pDONR223
NCBI SYMBOL	IL25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]