

CRISP3-cysteine-rich secretory protein 3 Gene

Reference: PTXBC069580

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Aeg2; CRISP-3; CRS3; SGP28; dJ442L6.3; cysteine-rich secretory protein 3; specific granule protein (28 kDa); cysteine rich secretory protein 3		

For research use only.

GENE DESCRIPTION	cysteine-rich secretory protein 3
GENE ID	10321
ORIGIN SPECIES	Human
ACCESSIONS	BC069580
VECTOR	pDONR223
NCBI SYMBOL	CRISP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacattatccagtgctgtgttctgctgtggtgctgttccatctttccagcaaatgaagataagatccgcttttactgctttgttaaccacccaaacacaaagtgcaggaggagattgtgaataagcacaatgaactgaggagagcagtatctccccctgccagaacatgctgaag
atggaatggaacaaaggctgcagcaaatgcccagaagtgggcaaacacagtgcaattacagacagtaacccaaaggatcgaatgacaagtctaaatgtgtgagaatctctacatgcaagtgcctccagctcatggtcacaagcaatccaaagctgtttgatgagtacaatgattt
tgactttgtgtaggccaaagactccaacgcagtggttgacattatacacaggtgtttgggtactcttcatacctcgttggatgtggaatgcctactgtcccaatcaaaagtctaaataactactatgtttgccaatattgtcctgctgtaattgggctaatagactatatgtccctatgaac
aaggagcacctgtgcccagttgccagataactgtgacgatggactatgcaccaatggtgcaagtacgaagatcttatagtaactgtaaaagttgaagctcacattaacctgtaaacatcagttggtcaggagacagttgcaaggcctcctgcaattgttcaaacagcatttattaa