



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MC3R-melanocortin 3 receptor Gene

Reference: PTXBC069599

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BMIQ9; MC3; MC3-R; OB20; OQTL; melanocortin receptor 3; obesity quantitative trait locus; melanocortin 3 receptor		

For research use only.

GENE DESCRIPTION	melanocortin 3 receptor
GENE ID	4159
ORIGIN SPECIES	Human
ACCESSIONS	BC069599
VECTOR	pDONR223
NCBI SYMBOL	MC3R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr