



Reference: PTXBC069601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LMX1.2; NPS1; LIM homeobox transcription factor 1-beta; LIM/homeobox protein 1.2; LMX-1.2; LIM homeobox transcription factor 1 beta		

For research use only.

GENE DESCRIPTION	LIM homeobox transcription factor 1, beta
GENE ID	4010
ORIGIN SPECIES	Human
ACCESSIONS	BC069601
VECTOR	pDONR223
NCBI SYMBOL	LMX1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]