



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

AQP10-aquaporin 10 Gene

Reference: PTXBC069607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AQPA_HUMAN; aquaporin-10; AQP-10; aquaglyceroporin-10; small intestine aquaporin; aquaporin 10		

For research use only.

GENE DESCRIPTION	aquaporin 10
GENE ID	89872
ORIGIN SPECIES	Human
ACCESSIONS	BC069607
VECTOR	pDONR223
NCBI SYMBOL	AQP10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctactcagcccccgcgtgaatacatggtgccacctccggatagcagcctctgcgcccgagctgctgcagagattctgggtgtgtttgtactcatgctctccaccaaggagctgtgccaccgctgtcaccagtggaagaacaaaggcaactctctaccatgttctgct
ggctctctgccgttagcatagaccatctacgtgggttgtaactgtcagggggccacctgaatcagcctctcctcgccatgtgcatgttgagacgctccccgggtcaagctccccattacatcttggtgcagttgtgtctgttctgtgttcgggagccacatgttctctaccatga
tgccctacagaacaatcaggtgggaacctgcagactgactggcccaagagacagcctcatttttgcacatcatctgccccctatctctcctgaacaattgctctccggatcaggttctgggacgcgagctgtattggggctctggccatctggacagacggaacaaggag
tccctgcgggtctggagcctgtgggtgtggtggatgctgacctgtgcctctgttaccgttgcggatccgggtgacactctctctaccatctgctgtgctgggtctgaagctctcagctgtgctgtaatgctgtgctgtaatgctgtgggtgggtgctgct

More info: Please contact us for any question: contact@proteogenix.fr