

C5orf20-chromosome 5 open reading frame 20 Gene

Reference: PTXBC069608

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C5orf20; dendritic cell nuclear protein 1; dendritic cell associated nuclear protein		
For research use only.			

GENE DESCRIPTION	chromosome 5 open reading frame 20
GENE ID	140947
ORIGIN SPECIES	Human
ACCESSIONS	BC069608
VECTOR	pDONR223
NCBI SYMBOL	C5orf20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcattacggagcagcaacccacatacagaactcgagaagccacggcctggagactgtacctggacacaaagactggagagaggagctgggtggggaacccagagttccagggtgccactccccagctccaccagagaacttgggaatgagctgctgccctgagtgccct
ctccaggccctcagtgagggtctctacccctccaggaggaacaaacctggccagctgggtgctcgagagggggcagttcaatctctccacaggggactctgcaactccaattctcagtgagcactctgcgagggccctcaggaccaggaatgaactgcatagcagcagaagga
agacagggccagaccagggcggaggagcccggaacatctggtttgtagttcagactctacccgttcacagttcacacagttctcaccgggaactcacacttgccctgtaccaagttttaaggcagtttaagctctgccatccgagacttcattttcttgatgaaaaatcactgaaatca
tcagatccatggcaccacacttcaattctcctaacagctggaaccgtcaggtggttcaggcctggtcttcgacttgatcactatccctcacctgtcagacagccagagcaggcgggtgagttcctcccaaacctccactgcattctcagctccaccgcaggcagccac
gtgcctgagtga