

HIST1H4A-histone cluster 1, H4a Gene

Reference: PTXBC069654

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4FA; histone H4; H4 histone family, member A; histone 1, H4a; histone cluster 1, H4a; histone cluster 1 H4 family member a		

For research use only.

GENE DESCRIPTION	histone cluster 1, H4a
GENE ID	8359
ORIGIN SPECIES	Human
ACCESSIONS	BC069654
VECTOR	pDONR223
NCBI SYMBOL	HIST1H4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggacgtgtaagggcgggaagggtttggtaaggggggtgccaagcgccaccgcaaggtgttcgtgacaacatccagggcatccaagccggccatccggcgtctgccccggcgtggcgggtgtgaagcggatctctgtctgatctacgaggagactcgggggtgct
caaggtgttttggagaacgtgatccgtgacgtgtcacctatacggagcagccaagcgcaagacagtcactgccatggacgtggtctacgcgttaagcgccaggagcacccttatggctttggcggttaa