



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GIP-gastric inhibitory polypeptide Gene

Reference: PTXBC069663

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

gastric inhibitory polypeptide; glucose-dependent insulinotropic polypeptide; incretin hormone

For research use only.

GENE DESCRIPTION	gastric inhibitory polypeptide
GENE ID	2695
ORIGIN SPECIES	Human
ACCESSIONS	BC069663
VECTOR	pDONR223
NCBI SYMBOL	GIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgggcccacgaagaccttctgctgctgctgctgctccctgttctctgctggcagtgaggactagagagaagaagagggtcacttcacgcgctctccccctcctgctctgttgatctcatgtctaaaggtagagagccctcaacctcaggccccaggtagcgggaaggcgctttcatcagtgacta
cagttatgcttcagtcagaagattcaccacaagaagcttgcagctctgctgctggcccaaaagggggaagaagaatgactgctggcgaacacaacatcaccanagagggtctcggggcgtgctggctggccactcaagtataatgaagaggagagcagtcaggccacagagctccca
ccgaagaaccccccagctataaagattgtctgcggcgaactctatcaacgaagctctgctcctctctgacacacaacactctcaggctcagctctctgtgta